

**Head Office / Manufacturer  
Elcam Medical**

BarAm 1386000,  
Israel  
Tel: 972 4 698 8120/1  
Fax: 972 4 698 0777  
sales@elcam.co.il

**USA  
Elcam Medical Inc.**

2 University Plaza, Suite 620  
Hackensack, NJ 07601  
Tel: +1 201 4571120  
Fax: +1 201 4571125  
info@elcam-medical.com

**Europe  
Elcam Medical Europe**


Via Emilia Romagna 15,  
41012 Carpi (Modena) Italy  
Tel: +39 059 697710  
Fax: +39 059 697705  
sales.eu@elcam-medical.com

**China  
Elcam China Office**

Room 3309, Jinhai Mansion  
No. 89 of Anshun Road,  
Shanghai, PRC. Zip: 200051  
Tel: +86 21 62829967  
Fax: +86 21 62829967  
johnliu02@126.com



 Elcam Devices  Elcam Stopcocks  Elcam Components



 **Elcam**  
**MEDICAL**  
Where everything connects

## Patient Monitoring



 [Elcam Devices](#)    [Elcam Stopcocks](#)    [Elcam Components](#)

# Patient Monitoring

Leverage your success to extend your competitive advantage with Elcam Medical's Patient Monitoring supreme solutions for this application.

Our unique patient monitoring products are the preferred solutions by industry leaders in this area. Their preference is based on our variety of solutions that truly meet theirs and their customers' needs.

Not just this, it is also based on unparalleled, end to end customer experience, tailored to fit your exact requirements in order to build your success in the marketplace.



## Sense-IT™ (DIPT)

**Disposable Integrated Pressure Transducer** Continuous Pressure Monitoring, is a mission-critical application.

Here you cannot allow your product, whether an Arterial Line, Central Venous Pressure line, or other pressure monitoring products, to fail. Elcam's Sense-IT™ (DIPT), is the first and only fully integrated pressure transducer. This integration provides you with the reliability and robustness no other transducer can provide, along with the space saving and ease of use so much needed in Intensive Care Units and Operating Rooms.

As in many of Elcam's products, Elcam's Sense-IT™ (DIPT) is also adaptable in dimensions such as tube fitments and cables connectors, to perfectly meet your requirements.

## DPT

**Disposable Pressure Transducer**

Sometimes you need just a straight transducer,

in applications where you do not need calibration or flushing, just plain pressure monitoring.

Elcam's DPT is simply the product for those tasks, no bells and whistles, but with Elcam's well known quality and service.



## Marvelous™ (MRVLS)

Specially designed for arterial and central venous blood sampling, this solution will give you the differentiating factor really needed for pressure monitoring lines.

This innovative product is named after the abbreviation - Minimal Residual Volume Luer Activated Swabbable Stopcock. It is a closed stopcock, integrated with a Luer Activated Valve, with a proprietary flow channel feature, allowing your customers to automatically flush blood residuals following blood sampling without the need to expose the system to a contaminated environment.



## Clear-IT™

### Closed Blood Sampling System

The Clear-IT™ is an advanced tool designed by Elcam in response to two common blood sampling needs: ICU Anemia due to blood waste associated with blood sampling, and infections due to exposure of patient's blood system to the atmosphere. The Clear-IT™ combines an In-Line Reservoir for waste collection and a stopcock with Marvelous™ proprietary handle that enables flushing the line and syringe tip with minimal residual volume.



## Flush Device

Sometimes you have an application where you need to maintain a very low flow.

It might be an application that requires keeping the vein open. Maybe you are producing your own pressure transducer and you are looking for a stand-alone flush device as a complimentary to your product line. Leveraged by Elcam's years of experience in pressure monitoring, Elcam's Flush Device is just the right solution for such applications. It allows you to be competitive and flexible, exactly where you need to be.