

Elcam Stopcocks

STOPCOCK WITH EXTENSION TUBE LINE

Elcam Medical offers its expertise in supplying stopcocks as cost effective sub-assemblies.



Your source of consistent quality and care



ELCAM STOPCOCK WITH EXTENSION TUBE LINE

As part of our continuous efforts to expand product range to best answer your stopcocks needs, Elcam is now offering its stopcocks with a variety of extension tube lines as cost effective sub-assemblies

Benefits

- Cost effective
- Compliance with highest quality standards
- Wide range – to meet market demands

Applications

- IV lines
- Introducers
- Central Venous Catheter
- Pressure Monitoring lines

Available configurations

Parameter	Values
Tube length (mm)	100-1500
Tube materials	PVC not made with DEHP products are available
Diameter O.D (mm)	2.8/4.1 (other ODs can be evaluated)
Connectors	Male Luer Locks (other connectors can be evaluated)
Covers	Optional

Tube fitment & Bonding

Elcam offers the widest variety of tube fitments for bonding, or oversizing, off the shelf or tailor made. Elcam offers its professional support in building and tweaking your bonding process.

Quality in every aspect:

Elcam's highly automated processes ensure consistent production of high quality products.

We operate in strict compliance with the international quality standards.

We bond so you won't have to!



www.elcam-medical.com

Head Office / Manufacturer Elcam Medical

BarAm 1386000,
Israel
Tel: 972 4 698 8120/1/2
Fax: 972 4 698 0777
sales@elcam.co.il

USA Elcam Medical Inc.

2 University Plaza, Suite 620
Hackensack, NJ 07601
Tel: +1 201 4571120
Fax: +1 201 4571125
info@elcam-medical.com

Europe Elcam Medical Europe

Via Emilia Romagna 15,
41012 Carpi (Modena) Italy
Tel: +39 059 697710
Fax: +39 059 697705
sales.eu@elcam-medical.com

China Elcam China Office

Room 3309, Jinhai Mansion
No. 89 of Anshun Road,
Sanghai, PRC. Zip: 200051
Tel: +86 21 62829967
Fax: +86 21 62829967
johnliu02@126.com