

Elcam Stopcocks

STOPCOCK WITH TUBE ASSEMBLY

Elcam Medical offers its expertise in supplying stopcocks as cost effective sub-assemblies.



Your source of consistent quality and care



STOPCOCK WITH TUBE ASSEMBLY

Elcam can bond your specified tubing to one of our stopcocks using either hand or automated processes. We utilize both solvent and UV bonding methods.

Benefits

- Cost effective
- Compliance with the highest quality standards
- Avoids the bonding process
- Flexibility – to bond your specified tubing

Applications

- IV lines
- Central Venous Catheter
- Pressure Monitoring lines

Available configurations

The available tubing lengths are between 1.5" and 25" with diameters between .110" (2.8mm) and .161" (4.1mm)

Tube fitment & Bonding

Elcam offers the widest variety of tube fitments for bonding, or oversizing, off the shelf or tailor made. Elcam offers its professional support in building and tweaking your bonding process.

Quality in every aspect:

Elcam's highly automated processes ensure production of consistently high quality products. We operate in strict compliance with the international quality standards.

Elcam's Stopcocks with Tube Assembly are CE approved for non-sterile products and 510(k) cleared.



www.elcam-medical.com

Head Office / Manufacturer Elcam Medical

BarAm 1386000,
Israel
Tel: 972 4 698 8120/1/2
Fax: 972 4 698 0777
sales@elcam.co.il

USA Elcam Medical Inc.

2 University Plaza, Suite 620
Hackensack, NJ 07601
Tel: +1 201 4571120
Fax: +1 201 4571125
info@elcam-medical.com

Europe Elcam Medical Europe

Via Emilia Romagna 15,
41012 Carpi (Modena) Italy
Tel: +39 059 697710
Fax: +39 059 697705
sales.eu@elcam-medical.com

China Elcam China Office

Room 3309, Jinhai Mansion
No. 89 of Anshun Road,
Sanghai, PRC. Zip: 200051
Tel: +86 21 62829967
Fax: +86 21 62829967
johnliu02@126.com