

# Safe2 Rotator™ Luer Connector


The connection you can rely on

Easy attach – Firm lock – 360° Rotational flexibility - Easy detach

Good and safe connection of medical tubing is critical to safe and effective therapy. Tubing disconnections can result in serious and irreversible patient injury or even death.

When it comes to stopcocks, apart from being securely connected, they must be convenient to use and perform their primary functions – direct fluids and provide access to the administration set.

With this in mind, Elcam Medical has designed the new Safe2 Rotator™.



## Secured Connections



 Elcam Devices

 Elcam Stopcocks

 Elcam Accessories

# Safe2 Rotator™ Luer Connector

The connection you can rely on

The Safe2 Rotator™ spinning lock feature provides ease of use and an extra layer of passive protection to your most critical fluid managing applications, without altering the traditional usability and common practice:

- Prevents accidental disconnections
- Facilitates convenient approach to injection and sampling ports (for example stopcock's side port and Y-site) during fluid administration
- Enables rotation of coupled luer lock connections and stopcocks
- Reduces tube kinks and twists in the tubing set (Kinks in the tubing can reduce the flow rates or cause blood hemolysis during transfusion or sampling)
- Prevents leakage due to over-tightening

Available configurations of the connector (stand-alone version):

Tube OD (mm)	Body Material	NUT Material	NUT color
6.8	PC	PP	○ ● ● ● ●
4.1	PC	PP	○ ● ● ● ● ●
4.1	ABS	PP	○ ● ● ● ● ●



Especially convenient in use with needleless connectors providing rotational flexibility of the distal part of the set.

**Elcam Medical's Safe2 Rotator™ – secured connection along the entire treatment**

**The Safe2 Rotator™ Luer Connector is available in both stand-alone and stopcock configurations. Stopcock configuration is CE approved for non-sterile products and 510(k) cleared**



[www.elcam-medical.com](http://www.elcam-medical.com)

Head Office/Manufacturer  
Elcam Medical

BarAm 1386000,  
Israel  
Tel: 972 4 698 8120  
Fax: 972 4 698 0777  
sales@elcam.co.il

USA  
Elcam Medical Inc.

2 University Plaza, Suite 620  
Hackensack, NJ 07601  
Tel: +1 201 4571120  
Fax: +1 201 4571125  
info@elcam-medical.com

Europe  
Elcam Medical Europe

Via Emilia Romagna 15,  
41012 Carpi (Modena) Italy  
Tel: +39 059 697710  
Fax: +39 059 697705  
sales.eu@elcam-medical.com

China  
Elcam China Office

Room 3309, Jinhai Mansion  
No. 89 of Anshun Road,  
Shanghai, PRC. Zip: 200051  
Tel: +86 21 62829967  
Fax: +86 21 62829967  
johnliu02@126.com