

Elcam Stopcocks

STANDARD MANIFOLDS

Monoblock Design

Elcam produces a wide range of manifolds that meet the demands of I.V. therapy including highflow systems, angioplasty, anesthesia, and hemodynamic pressure monitoring.



Your source of consistent quality and care














STANDARD MANIFOLDS

Monoblock Design

Elcam's Manifolds are available in 2, 3, 4 and 5 gang models in a variety of handle colors according to your needs. The right port can be female or male with a luer lock nut.

Our Manifolds can be supplied in materials suitable for ETO and Gamma irradiation, and/or lipid based admixtures upon request.

						
Handle Type		High Profile	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
		Low Profile	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
		GXL Handle	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
		T - Handle	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Right End Connections		Male Luer Lock	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
		Male Luer Slip	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
		Female Luer	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

Our Manifolds can be supplied with or without vented or sealed port covers.

Body Material:

Polycarbonate or Polysulfone.

Raw material available: either clear (for EtO sterilization only) or tinted blue (for gamma).

Quality in every aspect:

Elcam's highly automated processes ensure production of consistently high quality products. We operate in strict compliance with the international quality standards. Elcam's Standard Manifolds are 510k cleared and CE approved.

Colors:

Subject to Manifolds configuration.

