

Elcam Stopcocks

STANDARD STOPCOCKS

Standard Pressure, Standard Flow Path

Elcam Medical produces the widest range of stopcocks that meet the demands of hemodynamic pressure monitoring, IV therapy including highflow systems, angioplasty, anesthesia and other medical applications.



Your source of consistent quality and care














STANDARD STOPCOCKS

Standard Pressure, Standard Flow Path

Elcam's Standard Stopcocks are available in 1,2,3,4 way configurations with various handle styles to match different required applications.

Our Stopcocks can be supplied in materials suitable for ETO and Gamma irradiation, and/or lipid based admixtures upon request.



			1 Way (90°/360°)	3 Way (180°)	4 Way (360°)
Handle Type		High Profile		●	●
		Low Profile		●	●
		GXL Handle		●	●
		T - Handle			●
		2 Arm Low	●		
		2 Arm High	●		
End connections		Luer Lock			
		short Luer Lock		●	●
		Long Luer Lock		●	●
		Luer Slip	●	●	●
		Tube Pocket	●	●	●

Our Stopcocks can be supplied with or without vented or sealed port covers. At Elcam Stopcocks you can find the widest variety of "off the shelf" stopcocks tubing ports, for a variety of tubes, at all ports.

Body Material:

Polycarbonate, Polysulfone or Tritan.
Raw material available: either clear (for EtO sterilization only) or tinted blue (for gamma).

Colors:

Subject to stopcock configuration.



Quality in every aspect:

Elcam's state of the art equipment and automated processes ensure production of consistently high quality products. We operate in strict compliance with the international quality standards. Elcam's Standard Stopcocks are CE approved for non-sterile products and 510(k) cleared.



www.elcam-medical.com

Head Office/Manufacturer
Elcam Medical

BarAm 1386000,
Israel
Tel: 972 4 698 8120/1
Fax: 972 4 698 0777
menachem@elcam.co.il

USA
Elcam Medical Inc.

2 University Plaza, Suite 620
Hackensack, NJ 07601
Tel: +1 201 4571120
Fax: +1 201 4571125
info@elcam-medical.com

Europe
Elcam Medical Europe

Via Emilia Romagna 15,
41012 Carpi (Modena) Italy
Tel: +39 059 697710
Fax: +39 059 697705
sales.eu@elcam-medical.com

China
Elcam China Office

Room 3309, Jinhai Mansion
No. 89 of Anshun Road,
Shanghai, PRC. Zip: 200051
Tel: +86 21 62829967
Fax: +86 21 62829967
johnliu02@126.com