



**Intensive care & Anesthesia**



 **Elcam Devices**

 **Elcam Stopcocks**

 **Elcam Accessories**

# Intensive care & Anesthesia

As a long-time provider of solutions in this field, Elcam Medical is well aware of the importance of high quality treatment and reliable monitoring in critical care.

As a professional in the medical device industry dealing with Intensive care & Anesthesia, you need a committed partner that understands your needs and problems so you will have peace of mind when sourcing devices.

Elcam Medical is your source for the widest variety of Intensive Care & Anesthesia products including our renowned Elcam Stopcocks brand. We will partner with you to provide differentiating and innovative solutions along with the everyday devices you need, with special emphasis on maximizing patient and caregiver safety.

## Elcam Stopcocks for Intensive care & Anesthesia

Professionals procuring medical disposable stopcocks for treating critical patients know it is a critical device that needs to be sourced from experts. Elcam's broad range of stopcocks are designed to meet the requirements of Intensive care & Anesthesia including high flow systems, hemodynamic pressure monitoring and other applications. Elcam's stopcock range includes our specialty E-Safe™ line of stopcocks

**specially developed in order to provide our customers with solutions to minimize stopcock contaminations and enhance patient safety.**



### Closed Swabbable Stopcock

A Closed Stopcock, integrated with a Swabbable Luer Activated Valve, allowing for maintenance of a fully closed system through the entire use of the stopcock, as well as eliminating needle stick injuries when using a stopcock.



### Marvelous™ (MRVLS)

The innovative minimal residual volume luer activated swabbable stopcock, is a closed stopcock with a proprietary flow channel feature. The Marvelous™ allows your customers to continuously flush drugs and blood residuals, without the need to expose the system to the open air, and assuring minimal residual volume.

# Elcam Devices for Intensive care & Anesthesia

Patient Monitoring in Intensive care and Anesthesia is a mission critical application. Elcam's products developed for this application are the preferred solutions by industry leaders in the field, providing the enhanced safety required for fluid management during patient monitoring



## **Sense-IT™ (DIPT) Disposable Integrated Pressure Transducer**

Elcam's Sense-IT™ (DIPT), is the first and only fully integrated pressure transducer. This integration provides reliability and robustness as no other transducer on the market, along with the space saving and ease of use so much needed in Intensive Care Units and Operating Rooms.


Elcam Offers its DIPT components, the low-profile transducer and flush device, as standalone products, for specialized applications or for further set assembly.



## **Clear-IT™ Closed Blood Sampling System**

The Clear-IT™ is an advanced tool designed by Elcam in response to two common blood sampling needs: ICU Anemia due to blood waste associated with blood sampling, and infections due to exposure of patient's blood system to the atmosphere.

The Clear-IT™ combines an in-line reservoir for waste collection and a stopcock with Marvelous™ proprietary handle that enables flushing the line and syringe tip with minimal residual volume.



## **SafePort™ Manifold**

Unique function manifold with patented dual flow option innovative side-port valves.

Induction of anesthesia during operation is involved with enormous responsibility and requires multiple skills and focus. By integrating the SafePort™ Manifold into your solutions, you will simplify the anesthesiologists' work during procedures.

The SafePort™ Manifold is also available as a closed system with integrated Swabbable Luer Activated Valves, for minimizing leakages and exposure to contamination.



## **Drip Chambers**

Whenever there is a need for fluid or blood administration, we have the right drip chamber, for a wide range of intensive care and Anesthesia applications.

Elcam's drip chambers are available in different diameters and tube ports, for pediatric and adult applications, for transfusion and infusion needs.

## Elcam Accessories for Intensive care & Anesthesia

When sourcing medical devices for critical applications you need a full range of products from a source you can count on, someone who can support you in every situation.

We offer a broad range of accessories available for you from one source, whether it is an injection site, a connector, a fluid filter, or a clamp, at any size and any quantity.

As an Elcam Medical brand, we are well aware that delivering products is not enough. You need excellent service and support, from delivery time to design changes, from biocompatibility to optimizing the tube bonding processes. This is exactly the level of customer experience you will find when sourcing from Elcam Accessories™.



### Head Office / Manufacturer Elcam Medical

BarAm 1386000,  
Israel  
Tel: 972 4 698 8120/1  
Fax: 972 4 698 0777  
sales@elcam.co.il

### USA Elcam Medical Inc.

2 University Plaza, Suite 620  
Hackensack, NJ 07601  
Tel: +1 201 4571120  
Fax: +1 201 4571125  
info@elcam-medical.com

### Europe Elcam Medical Europe

Via Emilia Romagna 15,  
41012 Carpi (Modena) Italy  
Tel: +39 059 697710  
Fax: +39 059 697705  
sales.eu@elcam-medical.com

### China Elcam China Office

Room 3309, Jinhai Mansion  
No. 89 of Anshun Road,  
Shanghai, PRC. Zip: 200051  
Tel: +86 21 62829967  
Fax: +86 21 62829967  
johnliu02@126.com