


# Elcam Stopcocks

## STANDARD MANIFOLDS

### Monoblock Design

Elcam produces a wide range of manifolds that meet the demands of I.V. therapy including highflow systems, angioplasty, anesthesia, and hemodynamic pressure monitoring.



Your source of consistent quality and care



# STANDARD MANIFOLDS

## Monoblock Design

Elcam's Manifolds are available in 2, 3, 4 and 5 gang models in a variety of handle colors according to your needs. The right port can be female or male with a luer lock nut.

Our Manifolds can be supplied in materials suitable for ETO and Gamma irradiation, and/or lipid based admixtures upon request.

						
Handle Type		High Profile	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
		Low Profile	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
		GXL Handle	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
		T - Handle	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Right End Connections		Male Luer Lock	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
		Male Luer Slip	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
		Female Luer	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

Our Manifolds can be supplied with or without vented or sealed port covers.

### Body Material:

Polycarbonate or Polysulfone.

Raw material available: either clear (for EtO sterilization only) or tinted blue (for gamma).

### Quality in every aspect:

Elcam's highly automated processes ensure production of consistently high quality products. We operate in strict compliance with the international quality standards. Elcam's Standard Manifolds are 510k cleared and CE approved.

### Colors:

Subject to Manifolds configuration.



[www.elcam-medical.com](http://www.elcam-medical.com)

Head Office/Manufacturer  
Elcam Medical

BarAm 1386000,  
Israel  
Tel: 972 4 698 8120/1/2  
Fax: 972 4 698 0777  
menachem@elcam.co.il

USA  
Elcam Medical Inc.

1 University Plaza, Ste. 511  
Hackensack, NJ 07601  
Tel: +1 201 4571120  
Fax: +1 201 4571125  
info@elcam-medical.com

Europe  
Elcam Medical Europe

Via Emilia Romagna 15,  
41012 Carpi (Modena) Italy  
Tel: +39 059 697710  
Fax: +39 059 697705  
sales.eu@elcam-medical.com

China  
Elcam China Office

Room 3309, Jinhai Mansion  
No. 89 of Anshun Road,  
Shanghai, PRC. Zip: 200051  
Tel: +86 21 62829967  
Fax: +86 21 62829967  
johnliu02@126.com