


Elcam Stopcocks

E-Click

Detent indication stopcock

Providing tactile indication of handle position

The **E-Click** is intended to improve level of confidence for Medical teams - with regards to the handle position, by providing additional tactile orientation indication.




Your source of consistent quality and care



Introducing a new handle design for improved ergonomics and a better grip.

Elcam's new E-Click stopcock features a tactile indication handle with distinct handle stop over every 45° handle turn.



Benefits

- Tactile click every 45° turn
- Convenient and easy to use
- Enhances confidence in handle position
- New ergonomic handle design
- Lipid resistant


Applications

- IV Therapy
- Anesthesia
- Patient monitoring
- Other flow control applications

Available configurations

- T handle
- Handle colors- Blue, Red, White
- Luer lock

While focusing on better ergonomics and ease of use, the new stopcock is part of Elcam's road map plan to continue and develop advanced stopcock designs that better serve the needs of our customers, increase treatment safety while keeping handling and use easy and simple for the care givers.



Note: All closed position at 45°



www.elcam-medical.com

Head Office/Manufacturer
Elcam Medical

BarAm 1386000,
Israel
Tel: 972 4 698 8120/1
Fax: 972 4 698 0777
sales@elcam.co.il

USA
Elcam Medical Inc.

1 University Plaza, Ste. 511
Hackensack, NJ 07601
Tel: +1 201 4571120
Fax: +1 201 4571125
info@elcam-medical.com

Europe
Elcam Medical Europe

Via Emilia Romagna 15,
41012 Carpi (Modena) Italy
Tel: +39 059 697710
Fax: +39 059 697705
sales.eu@elcam-medical.com

China
Elcam China Office

Room 3309, Jinhai Mansion
No. 89 of Anshun Road,
Shanghai, PRC. Zip: 200051
Tel: +86 21 62829967
Fax: +86 21 62829967
johnliu02@126.com