

Flush Device



For maintaining a very low flow

Flush devices are used for maintaining slow continuous flow of solutions when KVO (Keep Vein Open) is required, thus to minimize the risk of blood clotting at the catheter tip, this is done by maintaining a very low flow in the line. The Flush device enables momentary use at high flow rate, for priming the IV line or for flushing it when so required.

Applications

- Pressure monitoring lines
- KVO (Keep Vein Open) post catheterization procedures

Features

Two Flushing options	• Pressing "Squeeze" (one handed) • Pulling the "PigTail"
Two Flow Rate configurations	• 3ml/hr • 30ml/hr (for use with infusion pump)
Flush Rate	• 125 ml/min
Biocompatible	• According to ISO 10993-1

Configurations

	Flow Rate	Configuration
	3 ml/hr	"T" Configuration FFM (One entry and two exits)
	3 ml/hr	"Inline" Configuration FM (one entry and one exit)
	30 ml/hr	"T" Configuration FFM (One entry and two exits)

Other configurations are optional and subject to quantities.

Elcam's Flush Device is manufactured to meet the highest regulatory demands, utilizing state of the art technology including automatic assembly lines and laser drilling.



www.elcam-medical.com

Head Office/Manufacturer Elcam Medical

BarAm 1386000,
Israel
Tel: 972 4 698 8120/1
Fax: 972 4 698 0777
sales@elcam.co.il

USA Elcam Medical Inc.

1 University Plaza, Ste. 511
Hackensack, NJ 07601
Tel: +1 201 4571120
Fax: +1 201 4571125
info@elcam-medical.com

Europe Elcam Medical Europe

Via Emilia Romagna 15,
41012 Carpi (Modena) Italy
Tel: +39 059 697710
Fax: +39 059 697705
sales.eu@elcam-medical.com

China Elcam China Office

Room 3309, Jinhai Mansion
No. 89 of Anshun Road,
Shanghai, PRC. Zip: 200051
Tel: +86 21 62829967
Fax: +86 21 62829967
johnliu02@126.com