

High Flow Stopcocks



High Flow Path Ways

Designed mainly to provide higher flow rates or viscous fluid applications where required. Major applications are within rapid infusion, anesthesia and blood related procedures.

High Flow Stopcocks

High Flow Path Ways

Elcam's High Flow Stopcocks are available in 1, 3 and 4 way configurations with various handle styles and luer locks to match different required applications. Our High Flow Stopcocks can be supplied in materials suitable for ETO and Gamma irradiation, and/or lipid based admixtures upon request.



			1 Way (90°/360°)	3 Way (180°)	4 Way (360°)
Handle Type		High Profile			
		Low Profile			
		GXL Handle			
		T - Handle			
		1 Arm HF			
		2 Arm High			
End connections		Long Luer Lock			
		short Luer Lock			
		Safe2 Rotator™			
		Luer Slip			
		Tube Pocket			

Our Stopcocks can be supplied with or without vented or sealed port covers.

At Elcam Stopcocks you can find the widest variety of "off the shelf" stopcocks tubing ports, for a variety of tubes, at all ports.

Body Material:

Polycarbonate or Polysulfone.

Raw material available: either clear (for EtO sterilization only) or tinted blue (for gamma).

Quality in every aspect:


Elcam's highly automated processes ensure production of consistently high quality products.

We operate in strict compliance with the international quality standards.

Elcam's High Flow Stopcocks are 510(k) cleared.

Colors:

Subject to stopcock configuration.



www.elcam-medical.com

Head Office/Manufacturer Elcam Medical

BarAm 1386000,
Israel
Tel: 972 4 698 8120/1
Fax: 972 4 698 0777
sales@elcam.co.il

USA Elcam Medical Inc.

1 University Plaza, Ste. 511
Hackensack, NJ 07601
Tel: +1 201 4571120
Fax: +1 201 4571125
info@elcam-medical.com

Europe Elcam Medical Europe

Via Emilia Romagna 15,
41012 Carpi (Modena) Italy
Tel: +39 059 697710
Fax: +39 059 697705
sales.eu@elcam-medical.com

China Elcam China Office

Room 3309, Jinhai Mansion
No. 89 of Anshun Road,
Shanghai, PRC. Zip: 200051
Tel: +86 21 62829967
Fax: +86 21 62829967
johnliu02@126.com