


Standard Manifolds



Monoblock Design


Elcam Medical produces the widest range of manifolds that meet I.V therapy including highflow systems, angioplasty, anesthesia and hemodynamic pressure monitoring.

Standard Manifolds

Monoblock Design

Elcam's Manifolds are available in 2, 3, 4 and 5 gang models in a variety of handle colors according to your needs. The right port can be female or male with a luer lock nut.

Our Manifolds can be supplied in materials suitable for ETO and Gamma irradiation, and/or lipid based admixtures upon request.



Handle Type		High Profile				
		Low Profile				
		GXL Handle				
		T - Handle				
Right End Connections		Male Luer Lock				
		Male Luer Slip				
		Female Luer				
		Safe2 Rotator™				

Our Manifolds can be supplied with or without vented or sealed port covers.

Body Material:

Polycarbonate or Polysulfone.

Raw material available: either clear (for EtO sterilization only) or tinted blue (for gamma).

Quality in every aspect:

Elcam's highly automated processes ensure production of consistently high quality products.

We operate in strict compliance with the international quality standards.

Elcam's Standard Manifolds are 510(k) cleared.

Colors:

Subject to Manifolds configuration.



www.elcam-medical.com

Head Office/Manufacturer Elcam Medical

BarAm 1386000,
Israel
Tel: 972 4 698 8120/1
Fax: 972 4 698 0777
sales@elcam.co.il

USA Elcam Medical Inc.

1 University Plaza, Ste. 511
Hackensack, NJ 07601
Tel: +1 201 4571120
Fax: +1 201 4571125
info@elcam-medical.com

Europe Elcam Medical Europe

Via Emilia Romagna 15,
41012 Carpi (Modena) Italy
Tel: +39 059 697710
Fax: +39 059 697705
sales.eu@elcam-medical.com

China Elcam China Office

Room 3309, Jinhai Mansion
No. 89 of Anshun Road,
Shanghai, PRC. Zip: 200051
Tel: +86 21 62829967
Fax: +86 21 62829967
johnliu02@126.com