
Elcam MEDICAL

YOUR SOURCE OF CONSISTENT QUALITY & CARE

Disposable Pressure Transducer (DPT)



Indicated for pressure monitoring, without calibration or flushing, because sometimes you need just a straight transducer

Patient Monitoring

Devices Stopcocks Accessories

Benefits

Ease of use

- Full visibility of the flow path
- Easy and safe priming and debubling
- Improved luer locks - utilizing Elcam's "ramp lock" patent for excellent connectivity

Customization

- Variety of configurations, colors, ports and connectors
- Optional – your logo printing

Reliability & Safety

- Elcam's High quality standards

Applications

For any pressure monitoring application where no flushing is required, such as:

- CathLab
- Intra Cranial Pressure
- Intra-abdominal pressure
- Urodynamic pressure
- Intrauterine pressure
- For further assembly

Available configurations

Female or Male Luer on either side, with or without Luer locks

Tube Sockets

Bonded Stopcocks
The DPT is available with any combination of above 3 options per customer request

Available Connectors:

Type 1	
Type 2	
Type 3	
Type 4	
Type 5	

Technical specifications:

Sensitivity	5µV/V/mmHg ±1%
Zero offset	Max offset 30 mmHg
Sterilization	EtO sterilization compatible
Shelf life	5 years from assembly date

Regulatory status

- Biocompatible according to ISO 10993-1
- Compliance with the following international standards: ANSI/AAMI BP22-1994 IEC 60601-1, IEC 60601-2-34

Intended use

- Available only as an OEM product
- Intended for single use only
- Maximum use duration is 96 hours

Elcam MEDICAL
Where everything connects

www.elcam-medical.com

Head Office/Manufacturer
Elcam Medical
BarAm 1386000,
Israel
Tel: 972 4 698 8120/1
Fax: 972 4 698 0777
sales@elcam.co.il

USA
Elcam Medical Inc.
1 University Plaza, Ste. 511
Hackensack, NJ 07601
Tel: +1 201 4571120
Fax: +1 201 4571125
info@elcam-medical.com

Europe
Elcam Medical Europe
Via Emilia Romagna 15,
41012 Carpi (Modena) Italy
Tel: +39 059 697710
Fax: +39 059 697705
sales.eu@elcam-medical.com

China
Elcam China Office
Room 3309, Jinhai Mansion
No. 89 of Anshun Road,
Shanghai, PRC. Zip: 200051
Tel: +86 21 62829967
Fax: +86 21 62829967
johnliu02@126.com